

Drive and stirring technology for

PHARMA | FOOD | COSMETICS | CHEMICALS | AUTOMATION




Laboratory agitator eco-Lab




- Steppless speed setting with scale 1-12
- Supply air connection $\varnothing 6\text{mm}$
- Exhaust air connection $\varnothing 10\text{mm}$
- Housing aluminium anodised
- ATEX certification
Ex II 2G Ex h T5 Gb
Ex II 2D Ex h T100°C Db
- Acetone-resistant seals
- Oil and silicone-free
- Protection class IP67
- Quick-action chuck for agitator shafts
 $\varnothing 1,5\text{mm} - \varnothing 13\text{mm}$
- Equipped with laboratory tripod with two cross clamps optional
- Equipped with drive shaft and impeller optional

WE GENERATE MOVEMENT IN SPECIAL ENVIRONMENTS.



DIMENSIONS :
LABORATORY TRIPOD :



| Variants | Art. No. | Rotation speed | Torque | Air consumption at 100 rpm and medium performance | A |
|---------------|----------|----------------|--------|---|---------|
| eco-Lab 450 | 5063 | 0-300 rpm | 4 Nm | 25 l/min | 50.5 mm |
| eco-Lab 450 G | 5073 | 0-600 rpm | 2 Nm | 25 l/min | 93 mm |

WE GENERATE MOVEMENT IN SPECIAL ENVIRONMENTS.