
Drive and stirring technology for

PHARMA | FOOD | COSMETICS | CHEMICALS | AUTOMATION




## Laboratory agitator eco-Lab INOX



WE GENERATE MOVEMENT IN SPECIAL ENVIRONMENTS.

## DIMENSIONS All dimensions in mm.



Variants	Art. no.	Speed	Torque	Air consumption at medium power	А
eco-Lab 450	5063	0-300 1/min	4 Nm / 35 lbf in	25 l/min	50.5 mm / 1.988 in
eco-Lab 450 G	5073	0-600 1/min	2 Nm / 18 lbf in	25 l/min	93 mm / 3.66 in
eco-Lab 900	5064	0-300 1/min	8 Nm / 71 lbf in	50 l/min	74.5 mm / 2.9 in
eco-Lab 900 G	5074	0-600 1/min	4 Nm / 35 lbf in	50 l/min	107 mm / 4.21



WE GENERATE MOVEMENT IN SPECIAL ENVIRONMENTS.